# **Electricité**

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# **Electricity**

*Three-phase sinusoidal mode*

Ref :  
**302 189**

Français – p 1

English – p 4

Version : 8009

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# 1 Asynchronous three-phase motor

## 1.1 Objective

- To measure the power consumed by an asynchronous three-phase motor in a star-delta connection.
- To measure a frequency of rotation.
- To calculate the sliding  $g = \frac{n_s - n}{n_s}$ .

## 1.2 Equipment

- An extra-low voltage three-phase generator;
- An extra-low voltage three-phase asynchronous motor;
- A double-pole change-over switch;
- A set of digital multimeters;
- An extra-low voltage wattmeter WJ2000;
- Connecting wires;
- A stroboscope.

## 1.3 Experiment

- ✓ Make the star-delta connection (for this refer to the simulation performed on the extra- low voltage three phase terminal board).
- ✓ Adjust the extra-low voltage three-phase generator:  $3 \times 9 \text{ V} - 50 \text{ Hz}$ .
- ✓ Observe the rotational movement of the rotor. Interchange the two phase wires; observe the change in direction of the rotation of the rotor.
- ✓ Measure for each connection:

- Line current  $I_Y$  or  $I_\Delta$
- The voltage
- The total power consumed  $P_Y$  or  $P_\Delta$
- The rotational speed of the rotor  $n_Y$  or  $n_\Delta$
- The sliding  $g_Y$  or  $g_\Delta$ .

For measuring the power, refer to the simulation performed on the extra-low voltage three-phase terminal board.

To measure the rotational speed  $n$  of the rotor using the stroboscope: increase the frequency of flashes; the rotational speed corresponds to the largest frequency for which the reference frame is immobile.

Check that the star circuit is balanced:  $I_N \approx 0$  A.

Check that the motor rotates slower in star connection.

Verify the relations:  $\frac{I_\Delta}{I_Y} = 3$  and  $\frac{P_\Delta}{P_Y} = 3$ .

Compare  $P$  measured with the value of  $P$  given by the relation:

$$P = U \times I \times \sqrt{3} \times \cos \varphi.$$

#### Examples of results:

$$n_S = \frac{f}{p} = \frac{53}{2} = 26,5 \text{ rps}$$

$$g = \frac{n_S - n}{n_S}$$

$n_S$  : rotational frequency of the rotating field

$f$  : Frequency of the current

$p$  : number of pole pairs

$g$  : sliding

$n_S$  : rotational frequency of the rotating field

$n$  : rotational frequency of the rotor

Connection	$I$ measured	$U$ measured	$P$ measured	$n$ measured	$g$
Star	78 mA	15.6 V	1.06 W	1378 rpm 23 rps	13%
Delta	171 mA	15.6 V	*	1554 rpm 25.9 rps	2%

\* The extra-low voltage wattmeter WJ2000 enables measurements only between the neutral and phase. (star connection).

$$\cos \varphi = \frac{P}{S} \qquad \cos \varphi = \frac{1,06}{15,6 \times 0,078 \times \sqrt{3}} \approx 0.50$$

## 2 Single phase asynchronous motor

### 2.1 Objective


To measure the power consumed by a single-phase asynchronous motor.

## 2.2 Equipment

- An extra-low voltage A.C. supply
- An extra-low voltage three-phase asynchronous motor
- A 33  $\mu\text{F}$  un-biased capacitor
- Two digital multimeters
- An extra-low voltage wattmeter WJ2000
- Connecting wires.

## 2.3 Experiment

Carry out the set-up as shown below

Two coils are set out of phase using a capacitor.

( Put on the switch; observe the direction of rotation of the rotor. Once the motor has started, remove the capacitor.

Measure the supply voltage, the line current and the power consumed; compare the power measured using the digital wattmeter to the value of the power given by the relation:

$$P = U \times I \times \cos \varphi .$$

## 3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**0825 563 563**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries